



# Innovative Inventory Management



## In Denmark, they know how to handle the cold...

Tenak is the most innovative supplier of cold storage inventory systems in the market. In addition to all the quality racks and boxes for freezers shown on the following pages, Tenak offer:

### Innovative racking systems:

- Cardboard racks
- Aluminium racks
- 'Side-up' racks to use your freezer's wasted space
- Adjustable-shelf racks and
- Custom-designed racks

### Innovative boxes:

- Origami flat-packed cryoboxes you can assemble in 30seconds Save space and save money
- PP cryoboxes for 15 and 50mL centrifuge tube storage

### Innovative tube racks:

- Origami tube racks for 10, 15 and 50mL centrifuge tubes

## ...and, it's easy to find the products you need:

Go to the comprehensive Product Reference guide at [www.tenak.com](http://www.tenak.com) and use the drop-down menu to easily identify the appropriate racks and boxes for your freezer, whatever the make/model or just ask us to do it for you.

Contact **CapellaScience** to maximise your storage capacity at economical prices.



Phone (East) 0401 888 088  
(West) 0423 533 245  
Email [enquiries@capellascience.com.au](mailto:enquiries@capellascience.com.au)



Biocontainment ★ Biodecontamination ★ Culturing

## Racks for Upright Freezers

Standard Racks	Price	Comfort Racks	Price	Box height	Box arrangement (H x D = Total)
TE23120	\$116	TE24220	\$195	50 mm	3x3 = 9
TE23122	\$119	TE24222	\$201	50 mm	4x3 = 12
TE23124	\$122	TE24224	\$215	50 mm	5x3 = 15
TE23130	\$119	TE24230	\$198	50 mm	3x4 = 12
TE23132	\$122	TE24232	\$203	50 mm	4x4 = 16
TE23134	\$125	TE24234	\$215	50 mm	5x4 = 20
TE23136	\$144	TE24236	\$227	50 mm	6x4 = 24
TE23138	\$150	TE24238	\$238	50 mm	7x4 = 28
TE23150	\$137	TE24250	\$215	50 mm	3x5 = 15
TE23152	\$141	TE24252	\$227	50 mm	4x5 = 20
TE23154	\$150	TE24254	\$238	50 mm	5x5 = 25
TE23156	\$156	TE24256	\$250	50 mm	6x5 = 30
TE23160	\$114	TE24260	\$195	75mm	3x2 = 6
TE23170	\$116	TE24262	\$195	75 mm	2x3 = 6
TE23172	\$119	TE24264	\$201	75 mm	3x3 = 9
TE23174	\$127	TE24266	\$212	75 mm	4x3 = 12
TE23180	\$119	TE24272	\$201	75 mm	2x4 = 8
TE23182	\$118	TE24274	\$212	75 mm	3x4 = 12
TE23184	\$132	TE24276	\$224	75 mm	4x4 = 16
TE23186	\$147	TE24278	\$235	75 mm	5x4 = 20
TE23190	\$132	TE24284	\$212	75 mm	2x5 = 10
TE23192	\$147	TE24286	\$224	75 mm	3x5 = 15
TE23194	\$153	TE24288	\$235	75 mm	4x5 = 20



Comfort

Standard

### Standard vs. Comfort

Standard racks feature traditional side-loading design.

Comfort racks feature sliding shelves, for ease of access and ergonomic benefits.

### Info about freezer racks

For full dimensions and more information about these racks visit [www.tenak.com](http://www.tenak.com)

## Racks for Chest Freezers

Classic Racks	Price	Comfort Racks	Price	For Boxes	No. of Boxes
TE21448	\$83	TE24500	\$90	50 mm	3
TE21449	\$86	TE24502	\$95	50 mm	5
TE21450	\$87	TE24504	\$98	50 mm	6
TE21451	\$92	TE24506	\$101	50 mm	7
TE21440	\$95	TE24508	\$106	50 mm	8
TE21441	\$96	TE24510	\$109	50 mm	9
TE21452	\$98	TE24512	\$112	50 mm	10
TE21454	\$99	TE24514	\$116	50 mm	11
TE21456	\$100	TE24516	\$121	50 mm	12
TE21478	\$103	TE24518	\$124	50 mm	13
TE21430	\$106	TE24520	\$127	50 mm	14
TE21432	\$112	TE24522	\$130	50 mm	15
TE21455	\$83	TE24530	\$90	75 mm	2
TE21457	\$83	TE24532	\$95	75 mm	3
TE21458	\$87	TE24534	\$98	75 mm	4
TE21460	\$90	TE24536	\$101	75 mm	5
TE21462	\$95	TE24538	\$106	75 mm	6
TE21464	\$101	TE24540	\$109	75 mm	7
TE21466	\$103	TE24542	\$112	75 mm	8
TE21467	\$109	TE24544	\$120	75 mm	9



Classic



Comfort

### Classic vs. Comfort

Classic racks feature a locking rod. Comfort racks feature Tenak's unique safety spring lock.

We have featured racks for 50mm and 75mm boxes: racks are also available for 32, 85, 100 and 130mm boxes.

## Racks for LN<sub>2</sub> Storage Full Size Boxes (Grid 9x9 or 10x10)

Classic Racks	Price	Comfort Racks	Price	For Boxes (Box height)	Number of Boxes
TE21448	\$83	TE24500	\$90	50 mm	3
TE21449	\$86	TE24502	\$95	50 mm	5
TE21450	\$87	TE24504	\$98	50 mm	6
TE21451	\$92	TE24506	\$101	50 mm	7
TE21440	\$95	TE24508	\$106	50 mm	8
TE21441	\$96	TE24510	\$109	50 mm	9
TE21452	\$98	TE24512	\$112	50 mm	10
TE21454	\$99	TE24514	\$116	50 mm	11
TE21456	\$100	TE24516	\$121	50 mm	12
TE21478	\$103	TE24518	\$124	50 mm	13
TE21430	\$106	TE24520	\$127	50 mm	14
TE21432	\$112	TE24522	\$130	50 mm	15



### Classic vs. Comfort

Standard racks feature a locking rod.

Comfort racks feature Tenak's unique Safety spring lock.

## Half Size Boxes (Grid 5x5)

Classic Racks	Price	Comfort Racks	Price	For Boxes (Box height)	Number of Boxes
TE21410	\$80	TE24490	\$99	50 mm	9
TE21412	\$85	TE24491	\$103	50 mm	10
TE21414	\$89	TE24492	\$108	50 mm	11
TE21416	\$92	TE24493	\$112	50 mm	12
TE21418	\$95	TE24494	\$115	50 mm	13
TE21420	\$98	TE24495	\$118	50 mm	14
TE21422	\$101	TE24496	\$121	50 mm	15



## Racks for Goblets

Codes	No. Goblets	Price
TE21922	2	\$112
TE21924	3	\$118
TE21926	4	\$124
TE21928	5	\$130



## Other types of racks available

**Microtitre Plate Racks:** Racks for all types of microtitre plates to suit all freezers, in both Comfort and Standard/Classic styles.



**Cardboard and Aluminium Racks:** Low-cost racks, available in Standard formats in a variety of sizes for chest and upright freezers.

**Side-Up Racks:** Innovative racks to maximise your storage capacity by filling the 'wasted' space in your freezer.



**Alu-flex Racks:** Aluminium racks with moveable shelves: versatile storage for all sizes of cryoboxes in chest freezers.

## Cardboard Cryoboxes, 133 x 133mm

Boxes/Dividers (height)	Box 50mm	Price/Pk 36	Box 75mm	Price/Pk 36
White box only	TE22005	\$89	TE22050	\$120
9x9 dividers only	TE22110	\$54	TE22134	\$87
10x10 dividers only	TE22120	\$54	TE22136	\$87
Boxes with dividers (height)	Box 50mm	Price/Pk 36	Box 75mm	Price/Pk 36
White, 9x9 dividers	TE22230	\$135	TE22235	\$205



\*All cardboard boxes also available in blue, red, yellow and green.

## Origami Boxes

Innovative 50mm flat-packed boxes save you space and money. Dividers available from 3x3 to 16x16 formats.

134 x 134mm footprint	Box 50mm	Price/Pk 100
Box, no dividers	TE22810	\$149
Box with 10x10 dividers	TE22820	\$249
Box with 9x9 dividers	TE22825	\$249
Box with 8x8 dividers	TE22830	\$249
Box with 7x7 dividers	TE22835	\$249

Scan to see how to assemble the Origami box in just 30 seconds or visit [www.capellascience.com.au/Origami](http://www.capellascience.com.au/Origami)



## Polypropylene Boxes

130 x 130mm, 9x9 grid	Box 50mm	Price/Pk 36
Clear	TE22310	\$169
Blue	TE22300	\$169
Green	TE22305	\$169
Red	TE22315	\$169
Yellow	TE22320	\$169



### Also Available

- Cardboard boxes in 32mm, 100mm and 130mm heights.
- Dividers from 3x3 to 16x16 grids
- Cardboard boxes with 136 x 136mm footprint
- Custom-printed boxes

## Polypropylene Boxes for Centrifuge Tubes

130 x 130 x 130mm	25x15mL tubes	Price/Pk 12	10x50mL + 2 x15mL	Price/Pk 12
Clear	TE22392	\$169	TE22362	\$169
Blue	TE22390	\$169	TE22360	\$169
Green	TE22391	\$169	TE22361	\$169
Red	TE22393	\$169	TE22363	\$169
Yellow	TE22394	\$169	TE22364	\$169



## Polycarbonate Cryoboxes

PC boxes for LN <sub>2</sub>	Code	Price/Pk
5x5 grid, 76x76x53mm, Red	TE22400	\$79 (Pack/8)
9x9 grid, 133x133x53mm, Assorted colours	TE22405	\$51 (Pack/4)
10x10 grid, 133x133x53mm, Blue	TE22410	\$51 (Pack/4)
9x9 grid, 133x133x96mm, Purple	TE22415	\$89 (Pack/4)



## Origami Tube Racks

Origami PP tube racks	Code	Price/Pack 5
For 45 x 12.5mm Ø tubes (9x5 array)	TE50510	\$49
For 28 x 18mm Ø tubes (7x4 array)	TE50530	\$49
For 15 x 30mm Ø tubes (5x3 array)	TE50550	\$49



Scan to see how to assemble the Origami box in just 30 seconds or visit [www.capellascience.com.au/Origami](http://www.capellascience.com.au/Origami)