

# THE SPACE RACE IS OVER YOU WIN!

## Increase your freezer capacity by up to 80%\*


Tenak's space-saving inventory gives you massive increases in your freezer storage capacity and drastic reductions in price per sample stored.

(\*by switching from 2.0mL traditional external-thread cryovials in 50mm cryoboxes to Tenak SlimTubes (0.5mL or 1.0mL) in Tenak 40mm cryoboxes and racks.)



### Tenak SlimTube Cryovials

Externally-threaded tubes which fit in a 10x10 array.



Product Description	Product Code
Tenak SlimTube 0.5mL Cryovials, Height 28.6mm, Starfoot, Labelled, Sterile, Pack/1,000 (10 bags x100)	TE78305**
Tenak SlimTube 1.0mL Cryovials, Height 36.1mm, Starfoot, Labelled, Sterile, Pack/1,000 (10 bags x100)	TE78310*
Tenak SlimTube 2.0mL Cryovials, Height 47.6 mm, Starfoot, Labelled, Sterile, Pack/1,000 (10 bags x100)	TE78315
Tenak SlimTube 4.0mL Cryovials, Height 73.3 mm, Starfoot, Labelled, Sterile, Pack /500 (5 bags x100)	TE78320
Tenak SlimTube 5.0mL Cryovials, Height 88.5mm, Starfoot, Labelled, Sterile, Pack /500 (5 bags x100)	TE78325

\*Fits in the Tenak 40mm cryoboxes

\*\*Fits in the Tenak 32mm and 40mm cryoboxes

## Tenak Cryoboxes

Storage Boxes 40mm size, grid 10x10



Product Description	Grid	Product Code
40mm White Cardboard Cryobox, Pack/36	None	TE22030
40mm White Cardboard Cryobox, Pack/36	9x9	TE22227
40mm White Cardboard Cryobox, Pack/36	10x10	TE22228

32mm, 50mm, 75mm, 100mm and 130mm cryoboxes also available.

## Tenak Comfort Freezer Racks

40mm size



HxD (no. of boxes)	Product Code	Style	Dimensions (mm, HxDxW)
4x4 = 16	TE24022	Comfort	171x562x140
5x4 = 20	TE24024	Comfort	214x562x140
6x4 = 24	TE24026	Comfort	256x562x140
7x4 = 28	TE24028	Comfort	299x562x140
8x4 = 32	TE24030	Comfort	341x562x140
9x4 = 36	TE24032	Comfort	384x562x140
4x5 = 20	TE24034	Comfort	171x683x140
5x5 = 25	TE24036	Comfort	214x683x140
6x5 = 30	TE24038	Comfort	256x683x140
7x5 = 35	TE24040	Comfort	299x683x140
8x5 = 40	TE24042	Comfort	341x683x140
9x5 = 45	TE24044	Comfort	384x683x140

Standard side-loading racks also available – ask for details