



Frontier™ 5000 Series

Centrifuge Bundles



Research Applications—Centrifuge Bundle Guide

Simplify your selection process with the OHAUS guide to centrifuge and rotor bundles! Featuring combinations designed around specific centrifugation applications, from high speed to cell culture to microplate, OHAUS has the right centrifuge and rotor packages to help you *Do More* in your laboratory.

**MASSIVE
SAVINGS**
on CENTRIFUGE
BUNDLES

Microcentrifuge Bundles

Designed for many biological applications to spin small liquid samples at high speeds. Our bundles include options for high throughput applications to biocontainment requirements.

Mid-Speed Micro - Non-Refrigerated Bundles



FC5513

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor	List Price (excl. GST)	Bundle Price (excl. GST)	
Standard	24	1.5/2.0 mL	13,500	17,317	Centrifuge FC5513	30370691		\$3,190	\$3,100	View brochure
					Rotor 24 x 1.5/2.0 mL					
Biocontainment	24	1.5/2.0 mL	13,500	17,317	Centrifuge FC5513	30393187		\$3,451	\$3,385	View brochure
					Sealable Rotor 24 x 1.5/2.0 mL	30472302				

High-Speed Micro - Non-Refrigerated Bundles



FC5515

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor	List Price (excl. GST)	Bundle Price (excl. GST)	
Standard	24	1.5/2.0 mL	15,200	21,953	Centrifuge FC5515	30130866		\$4,937	\$4,695	View brochure
					Rotor 24 x 1.5/2.0 mL	30130870				
Biocontainment	24	1.5/2.0 mL	15,200	21,953	Centrifuge FC5515	30130866		\$5,225	\$4,990	View brochure
					Sealable Rotor 24 x 1.5/2.0 mL	30130871				
Biocontainment	30	1.5/2.0 mL	13,500	19,150	Centrifuge FC5515	30130866		\$5,430	\$5,275	View brochure
					Sealable Rotor 30 x 1.5/2.0 mL	30130872				
High Volume	44	1.5/2.0 mL	13,500	17,317	Centrifuge FC5515	30130866		\$5,304	\$4,995	View brochure
					Rotor 44 x 1.5/2.0 mL	30130879				

High-Speed Micro - Refrigerated Bundles



FC5515R

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor	List Price (excl. GST)	Bundle Price (excl. GST)	
Standard	24	1.5/2.0 mL	15,200	21,953	Centrifuge FC5515R	30130868		\$9,207	\$8,500	View brochure
					Rotor 24 x 1.5/2.0 mL	30130870				
Biocontainment	24	1.5/2.0 mL	15,200	21,953	Centrifuge FC5515	30130868		\$9,495	\$8,800	View brochure
					Sealable Rotor 24 x 1.5/2.0 mL	30130871				
Biocontainment	30	1.5/2.0 mL	13,500	19,150	Centrifuge FC5515R	30130868		\$9,700	\$9,100	View brochure
					Sealable Rotor 30 x 1.5/2.0 mL	30130872				
High Volume	44	1.5/2.0 mL	13,500	17,317	Centrifuge FC5515R	30130868		\$9,574	\$8,780	View brochure
					Rotor 44 x 1.5/2.0 mL	30130879				

Cell Culture Bundles

Designed for many cell culture applications to separate or concentrate cells from culture media and/or to collect extracellular products. Our bundles include accessories for 15ml and/or 50ml conical tubes, ideal for your cell culture needs.

Low Volume – Cell Culture Bundles



FC5707

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	8 x 15ml RB/ 4 x 15ml FA	15mL	6,800	4,445	Centrifuge FC5707	30393189		\$2,320	\$2,270	View brochure
					Rotor Angle 8 x15ml					

Low Volume – Cell Culture Bundles



FC5706

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	12	15mL RB/FA	6,000	4,427	Centrifuge FC5706	30332131		\$2,995	\$2,750	View brochure
					Rotor Angle 12 x15ml					
Non-Refrigerated	6	15 or 50mL	6,000	-	Centrifuge FC5706	30130875		\$4,556	\$4,200	View brochure
				4,427	Rotor Angle 6 x 50 mL	30130878				
				4,225	Adapters (qty 6) 15 mL Tubes	30130893				

Mid Volume – Cell Culture Bundles



FC5714

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	8 x 50mL OR 12 x 15mL	50mL FA OR 15mL FA	4,000	-	Centrifuge FC5714	30314810		\$8,478	\$8,290	View brochure
				-	Rotor Swing Out 4 x100mL Sealable	30314823				
				2,486	Adapter Bucket, 3 x 15 mL (QTY. 4)	30314882				
				2,486	Adapter Bucket, 2 x 50mL (QTY. 4)	30314881				

Mid Volume – Cell Culture Bundles



FC5718



FC5718R

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	8 x 50mL OR 12 x 15mL	50mL FA OR 15mL FA	5,000	-	Centrifuge FC5718	30314812		\$10,413	\$10,090	View brochure
				-	Rotor Swing Out 4 x 100mL Sealable	30314823				
				3,885	Adapter Bucket, 3 x 15 mL (QTY. 4)	30314882				
				3,885	Adapter Bucket, 2 x 50mL (QTY. 4)	30314881				
Refrigerated	8 x 50mL OR 12 x 15mL	50mL FA OR 15mL FA	5,000	-	Centrifuge FC5718R	30314814		\$15,543	\$12,990	View brochure
				-	Rotor Swing Out 4 x 100mL Sealable	30314823				
				3,885	Adapter Bucket, 3 x 15 mL (QTY. 4)	30314882				
				3,885	Adapter Bucket, 2 x 50mL (QTY. 4)	30314881				

High Volume – Cell Culture Bundles



FC5816



FC5816R

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	28 x 15 mL OR 12 x 50 mL	50mL FA OR 15mL FA	4,500	-	Centrifuge FC5816	30314816		\$12,835	\$11,490	View brochure
				-	Rotor Swing Out 4 x 250mL	30314820				
				3,712	Adapter Rack, 7 x 15 mL (QTY. 4)	30314913				
				3,780	Adapter Bucket, 3 x 50mL (QTY. 4)	30314906				
Refrigerated	28 x 15 mL OR 12 x 50 mL	50mL FA OR 15mL FA	4,500	-	Centrifuge FC5816R	30314818		\$18,435	\$14,990	View brochure
				-	Rotor Swing Out 4 x 250mL	30314820				
				3,712	Adapter Rack, 7 x 15 mL (QTY. 4)	30314913				
				3,780	Adapter Bucket, 3 x 50mL (QTY. 4)	30314906				

High Volume – Cell Culture Bundles



FC5916



FC5916R

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	68 x 15 mL OR 28 x 50 mL	50mL FA OR 15mL FA	4,000	3,434	Centrifuge FC5916	30553036		\$17,222	\$16,990	View brochure
					Rotor Swing Out 4 x 750mL with buckets	30553085 & 30553121				
					Adapter Rack, 17 x 15 mL (QTY. 4)	30553134				
					Adapter Bucket, 7 x 50mL (QTY. 4)	30553130				
Refrigerated	68 x 15 mL OR 28 x 50 mL	50mL FA OR 15mL FA	4,500	4,346	Centrifuge FC5916R	30553101		\$22,382	\$21,790	View brochure
					Rotor Swing Out 4 x 750mL with buckets	30553085 & 30553121				
					Adapter Rack, 17 x 15 mL (QTY. 4)	30553134				
					Adapter Bucket, 7 x 50mL (QTY. 4)	30553130				

Microtitre Plate Bundles

Microtitre plates are commonly used for analytical research and clinical diagnostic testing. Our bundles include accessories for microtitre plates with both refrigerated and non-refrigerated options.



FC5718



FC5718R

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	6 x Microtiter plates	Microtiter plates	4,500	2,716	Centrifuge FC5718	30314812		\$9,625	\$9,390	View brochure
					Rotor Swing Out 2 x 3MTP	30314824				
Refrigerated	6 x Microtiter plates	Microtiter plates	4,500	2,716	Centrifuge FC5718R	30314814		\$14,755	\$12,600	View brochure
					Rotor Swing Out 2 x 3MTP	30314824				

Clinical Bundles

Ideal for the separation of whole blood components, such as serum, plasma, red blood cells, as well as other bodily fluids. Our clinical centrifuges can accommodate common blood draw tubes. Please see specifications below for more details.

Low Volume - Clinical Bundles



FC5707

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	8 x 15ml RB OR 4 x 15ml FA OR 8 x 5ml Blood Tubes OR 8 x 7ml Blood Tubes	15mL OR 5ml Blood Tubes (13mm x 75mm) OR 7ml Blood Tubes (13mm x 100mm)	6,800	-	Centrifuge FC5707	30393189		\$3,076	\$2,960	View brochure
				4,445	Rotor Angle 8 x15ml	30472305				
				3,618	Adapters (qty 8) 5ml Blood Tubes	30130890				
				4,290	Adapters (qty 8) 7ml Blood Tubes	30130889				

Low Volume - Clinical Bundles



FC5706

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	6 x 5ml Blood Tubes	5ml Blood Tubes (13mm x 75mm)	4,000	1,878	Centrifuge FC5706	30130875		\$3,573	\$3,190	View brochure
					Rotor Swing Out 6 x 5ml *Comes with buckets	30130880				

Mid Volume - Clinical Bundles



FC5714

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	16 x 5-7ml Blood Tubes OR 16 x 10ml Blood Tubes	5-7ml Blood Tubes (13mm x 75/100mm) OR 10ml Blood Tubes (16mm x 100mm)	4,000	2,418	Centrifuge FC5714	30314810		\$8,500	\$8,190	View brochure
					Rotor Swing Out 4 x100mL Sealable	30314823				
					Bucket with screw cap	30314861				
					Adapters Rack (qty 4) 5-7ml Blood Tubes	30314879				
					Adapters Rack (qty 4) 10ml Blood Tubes	30314880				

Mid Volume - Clinical Bundles



FC5718R

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Refrigerated	16 x 5-7ml Blood Tubes OR 16 x 10ml Blood Tubes	5-7ml Blood Tubes (13mm x 75/100mm) OR 10ml Blood Tubes (16mm x 100mm)	4,000	2,418	Centrifuge FC5718R	30314814		\$14,294	\$12,500	View brochure
					Rotor Swing Out 4 x 100mL Sealable	30314823				
					Bucket with screw cap	30314861				
					Adapters Rack 5-7ml Blood Tubes (QTY. 4)	30314879				
					Adapters Rack 10ml Blood Tubes (QTY. 4)	30314880				

High Volume - Clinical Bundles



FC5816

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	40 x 5-7ml Blood Tubes OR 32 x 10ml Blood Tubes	5-7ml Blood Tubes (13mm x 75/100mm) OR 10ml Blood Tubes (16mm x 100mm)	4,500	3,599 (5-7ml tube) 3,689 (4-10ml tube)	Centrifuge FC5816	30314816		\$13,155	\$12,500	View brochure
					Rotor Swing Out 4 x 250mL	30314820				
					Screw Cap (QTY. 4), O-rings	30314901 30314902				
					Adapter Bucket, 10 x 5-7mL (QTY. 4)	30314917				
					Adapter Rack, 8 x 4-10 mL (QTY. 4)	30304367				



FC5816R

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Refrigerated	40 x 5-7ml Blood Tubes OR 32 x 10ml Blood Tubes	5-7ml Blood Tubes (13mm x 75/100mm) OR 10ml Blood Tubes (16mm x 100mm)	4,500	3,599 (5-7ml tube) 3,689 (4-10ml tube)	Centrifuge FC5816R	30314818		\$15,362	\$14,700	View brochure
					Rotor Swing Out 4 x 250mL	30314820				
					Screw Cap (QTY. 4), O-rings	30314901 30314902				
					Adapter Bucket, 10 x 5-7mL (QTY. 4)	30314917				
					Adapter Rack, 8 x 4-10 mL (QTY. 4)	30304367				

High Volume – Clinical Bundles



FC5916

Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor/Adapters	List Price (excl. GST)	Bundle Price (excl. GST)	
Non-Refrigerated	108 x 5-7ml Blood Tubes OR 72 x 10ml Blood Tubes	5-7ml Blood Tubes (13mm x 75/100mm) OR 10ml Blood Tubes (16mm x 100mm)	4,000	3,434	Centrifuge FC5916	30553036		\$17,745	\$17,200	View brochure
					Rotor Swing Out 4 x 750mL with buckets & lids/O-rings	30553085, 30553121, 30553122 & 30553123				
					Adapters Rack 5-7ml Blood Tubes (QTY. 4)	30553139				
					Adapters Rack 10ml Blood Tubes (QTY. 4)	30553140				
Refrigerated	108 x 5-7ml Blood Tubes OR 72 x 10ml Blood Tubes	5-7ml Blood Tubes (13mm x 75/100mm) OR 10ml Blood Tubes (16mm x 100mm)	4,500	4,346	Centrifuge FC5916R	30553101		\$22,905	\$22,400	View brochure
					Rotor Swing Out 4 x 750mL with buckets & lids/O-rings	30553085, 30553121, 30553122 & 30553123				
					Adapters Rack 5-7ml Blood Tubes (QTY. 4)	30553139				
					Adapters Rack 10ml Blood Tubes (QTY. 4)	30553140				



FC5916R

Hematocrit Bundles

Designed to determine hematocrit values (HCT) for testing in biochemistry, immunity, genetics, blood separation and for general clinical tests.



Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor	List Price (excl. GST)	Bundle Price (excl. GST)	
Standard	24	Hematocrit Capillaries	13,000	16,058	Centrifuge FC5513	30393187		\$3,710	\$3,660	View brochure
					Rotor 24 x Capillaries	30472304				

PCR Bundle

Perfect for spinning PCR strips for PCR reaction set-ups.



Type	Max. Tube Qty	Volume	RPM	Max. RCF	Bundle Contents	Item No	Rotor	List Price (excl. GST)	Bundle Price (excl. GST)	
Standard	4x8-Place PCR Strips	0.2 mL PCR 8-strips	15,000	15,342	Centrifuge FC5515	30130866		\$5,724	\$5,650	View brochure
					Rotor 4x8-Place PCR Strips	30130874				