



VCE-800

Pharmacy



A Low Cost, High Quality Solution for Pharmacy Compounding

The Germfree VCE-800 is a Containment Ventilated Enclosure (CVE) as defined in USP <795> and USP <800>. The enclosure provides a work area for the safe handling of dust, powders and particulates associated with compounding. The blower provides between 50-70 cfm of inward airflow directing particulate toward the HEPA filter at the top of the unit. All exhausted air is HEPA filtered removing particulate at an efficiency of 99.99% for particulate 0.3 microns in size and returned to the room.



Equipment pictured may include custom or optional features

Features & Options

- HEPA Filter
 - Externally Placed to Reduce Vibration
- Can be Easily Leak Tested
- Configured for Easy Certification



Specifications

Standard features

- Meets USP <800> for Non-Sterile Hazardous Compounding
- Meets USP <795> for Non-Sterile Non-Hazardous Compounding
- Powder Coated Construction
Outer cabinet is constructed from powder coated steel and framed acrylic. Work surface is all powder coated steel
- HEPA Filter
High capacity motor blower system with speed control to extend HEPA filter life
- Ergonomics
Front viewing panel is top hinged and self-supporting for easy access to the work area for cleaning and equipment ingress

- Designed for independent certification

- Includes exhaust piece for external and redundant venting options to meet USP <800>

Electrical

- LED Lighting
- Separate lighted power ON/OFF indicator switches for blower and lighting
- Ten foot power cord with molded grounded plug
- Differential pressure gauge located in control panel

Model	Overall Dimensions		
	W	D	H
VCE-800-4	48"	23.188"	46.25"
VCE-800-3	36"	23.188"	46.25"

OPTIONAL FEATURES

- VCE-4-A8A36
4' Base Stand – 36" H
- VCE-A11
Casters for Base Stand
- RHF-4
Redundant HEPA Filter
- SBC
Stainless Steel Bin Cart
(30" W x 26-1/2" D x 29" H)
- BINS
Plastic Bins 14" x 6" x 6"
(Qty 12)

*Equipment pictured may include custom or optional features



Technical drawings

Please consult with a Germfree Sales Engineer for the most up to date engineering drawings.

(Page Intentionally Left Blank)

