



Kuhner shaker


# SCALE UP WITH KUHNER SB10-X NOW!

## KUHNER SB10-X SINGLE-USE BIOREACTOR

FREE 4 x PK3 Basic Bags  
with any SB10-X Purchase

Order an IQ/OQ with SB10-X  
purchase & receive an additional  
FREE PK3 Basic Bags

Order a SB10-X requalification  
& receive a FREE PK3 Basic Bags



Item Code	Description	UOM
106857	SB10-X	EA
105328	3 x basic single-use bags for SB10-X with no optical sensors	PK3
105332	3 x standard 3L single-use bags for SB10-X with optical chemosensors for pH/DO and tubes	PK3
105327	3 x standard single-use bags for SB10-X with optical chemosensors for pH/DO and tubes	PK3
105329	3 x perfusion single-use bags for SB10-X with 2 exhaust filters, optical chemosensors for pH/DO & additional tubes	PK3
107450	IQ/OQ OSB (SB10/50/200-X)	EA

### Kuhner SB10-X Orbital Shaken Benchtop Bioreactor for Single-Use Bags

- Cultivation of human, mammalian & plant cells in a single-use bag
- 4 different single-use bag options. Shaking speed 20rpm - 140rpm
- Temperature, pH & DO control, feeding & low shear stress
- Very few bubbles (no or less pluronic) & very little foam (no or less anti-foam)
- Similar hydrodynamics (between flasks & SB10X) & effective CO2 stripping

### Single-Use Bags

- Gamma-irradiated. Working volume from 4L up to 12L
- Quick set up times. No additional mixing device required
- Feeding, inoculation, harvesting & sampling ports incorporated



Kuhner SB10-X Brochure



Single-use bags for  
Kuhner OSB Brochure



Kuhner Notes



Orbital Shaken  
Bioreactors Brochure



☎ 02 9575 7512  
✉ sales@capellascience.com.au  
💻 capellascience.com.au

DOWNLOAD  
PROMOTION HERE



Contact us to learn more about this exclusive promotion on the Kuhner SB10-X Orbital Shaken Benchtop Bioreactor for Single-Use Bags. Promotion valid until 31.03.25. Does not apply to existing, pending, prior orders. Research use only and not for use in diagnostic procedures. Excludes GST & Freight. Images displayed are for visual purposes & may vary from the actual product. Call Capella Science now 02 9575 7512 to order or email sales@capellascience.com.au