

# STERIDIUM

## 160L Benchtop Oven

Model ob160



# 160L Benchtop Oven

## Model ob160

### Product Specifications

Chamber Capacity	160L
Chamber Dimensions	515 x 520 x 586 mm (w d h)
External Dimensions	658 x 727 x 857 mm (w d h)
Chamber Material	Stainless Steel (304 L)
External Material	Durable epoxy powder-coated steel
Temperature Range	(ambient+7) to 200C
Temperature Sensor	pt100
Temporal Stability	+/- 0.1C
Spatial Stability	+/- 2C
Shelf Positions	5
Shelves Supplied	2
Shelf Support Spacing	89 mm
Insulation Type	Mineral wool
Insulation Thickness	50 mm



**Brisbane, Australia**

+61 7 55 931 352

sales@steridium.com

steridium.com



# 160L Benchtop Oven

Model ob160

## Product Specifications

Electrical Supply	Single phase, 220-240V, 50Hz (alternatives available on request)
Electrical Connection Rating	10 Amps
Circulating fans	1 - Rear mounted
Countdown Timer	Included
Controller Type	Intelligent PID - microprocessor
Readout	Digital LCD
Place of manufacture	Brisbane, Australia



STERIDIUM

# 160L Benchtop Oven

## Model ob160

### Options / Accessories

- Additional Shelf
- Side Access Port - 50 mm
- Calibration Certificate
- 8 Step Programmable Controller - 4 Ramp/Dwell Pairs
- ISPM15 Compliant Export Crating
- Drying of moisture laden content - see model dg160
- Soil Drying - see model ds160



**Brisbane, Australia**

+61 7 55 931 352

[sales@steridium.com](mailto:sales@steridium.com)

[steridium.com](http://steridium.com)

