

STERID IUM

2400L General Drying Oven

Model dg2400



2400L General Drying Oven

Model dg2400

Product Specifications

Chamber Capacity	2400L
Chamber Dimensions	865 x 910 x 1434 mm (w d h)
External Dimensions	2220 x 1154 x 2110 mm (w d h)
Chamber Material	Stainless Steel (304 L)
External Material	Durable epoxy powder-coated steel
Temperature Range	(ambient+7) to 150C
Temperature Sensor	pt100
Temporal Stability	+/- 0.1 C
Spatial Stability	+/- 2C
Shelf Positions	30 (15 per side)
Shelves Supplied	8 (4 per side)
Shelf Support Spacing	76
Insulation Type	Mineral wool
Insulation Thickness	50 mm



Brisbane, Australia
 +61 7 55 931 352
 sales@steridium.com
 steridium.com



2400L General Drying Oven

Model dg2400

Product Specifications

Electrical Supply	Three phase, 380-41V, 50Hz (alternatives available on request)
Electrical Connection Rating	20 Amps / Phase
Circulating fans	2 - Base mounted
Countdown Timer	Included
Controller Type	Intelligent PID - microprocessor
Readout	Digital LCD
Place of manufacture	Brisbane, Australia
Moisture Extraction System	High performance drying of moisture-laden content utilising fresh air intake and hot air discharge
Exhaust	2 X 100 mm - Top Mounted




2400L General Drying Oven

Model dg2400

Options / Accessories

- Additional Shelf
- Modulating Valve
Increased maximum temperature to 200C
- Castors
- Side Access Port - 50 mm
- Calibration Certificate
- 8 Step Programmable Controller - 4 Ramp/Dwell Pairs
- ISPM15 Compliant Export Crating



Brisbane, Australia

+61 7 55 931 352

sales@steridium.com

steridium.com