

STERIDIUM

## 810L General Drying Oven

Model dg810



# 810L General Drying Oven

## Model dg810

### Product Specifications

|                       |                                   |
|-----------------------|-----------------------------------|
| Chamber Capacity      | 810L                              |
| Chamber Dimensions    | 842 x 705 x 1373 mm (w d h)       |
| External Dimensions   | 1048 x 881 x 1982 mm (w d h)      |
| Chamber Material      | Stainless Steel (304 L)           |
| External Material     | Durable epoxy powder-coated steel |
| Temperature Range     | (ambient+7) to 150C               |
| Temperature Sensor    | pt100                             |
| Temporal Stability    | +/- 0.1 C                         |
| Spatial Stability     | +/- 2C                            |
| Shelf Positions       | 16                                |
| Shelves Supplied      | 3                                 |
| Shelf Support Spacing | 77                                |
| Insulation Type       | Mineral wool                      |
| Insulation Thickness  | 50 mm                             |



# 810L General Drying Oven

## Model dg810

### Product Specifications

|                              |  |
|------------------------------|--|
| Electrical Supply            | Three phase, 380-41V, 50Hz<br>(alternatives available on request)                                  |
| Electrical Connection Rating | 16 Amps / Phase  |
| Circulating fans             | 1 - Base mounted   |
| Countdown Timer              | Included   |
| Controller Type              | Intelligent PID - microprocessor   |
| Readout                      | Digital LCD  |
| Place of manufacture         | Brisbane, Australia  |
| Moisture Extraction System   | High performance drying of moisture-laden content utilising fresh air intake and hot air discharge |
| Exhaust                      | 100 mm - Top Mounted   |




# 810L General Drying Oven

## Model dg810

### Options / Accessories

- Additional Shelf
- Modulating Valve  
Increased maximum temperature to 200C
- Castors
- Side Access Port - 50 mm
- Calibration Certificate
- 8 Step Programmable Controller - 4 Ramp/Dwell Pairs
- ISPM15 Compliant Export Crating



**Brisbane, Australia**

+61 7 55 931 352

[sales@steridium.com](mailto:sales@steridium.com)

[steridium.com](http://steridium.com)