



25% OFF OHAUS GENIE™ SOLUTIONS

THE POWERFUL COMBINATION OF OHAUS
AND GENIE™ - STRONGER TOGETHER



VORTEX-GENIE™

Code	Model	Control	Speed Range	Timer	Orbit	Mode	Display
83041399	Vortex-Genie™ 2 (SI-0297AU)	Analog	600–3200rpm	—	4mm	Touch/Hands Free	—
83041246	Vortex-Genie™ 2T (SI-T296AU)	Analog	600–3200rpm	1–60 sec	4mm	Touch/Hands Free	—
83041316	Digital Vortex-Genie™ 2 (SI-A296AU)	Digital	500–3000rpm	1–99 sec	3.5mm	Touch/Hands Free	Red LED
83041309	ABI Digital Vortex-Genie™ (SI-A596AU)	Digital	500–3000rpm	1–99 sec	3.5mm	Touch/Hands Free	Red LED
83041255	Vortex-Genie™ Pulse (SI-P296AU)	Digital	500–3000rpm	1–99 min	3mm	Touch/Hands Free	Red LED



MULTI-MODE ROTATOR GENIE™

Code	Model	Control	Rotating Speed	Rocking Speed	Tilt Angle	Motion	Capacity
83041253	Roto-Rock Genie™ (SI-R106AU)	Touch	1–35 rpm	1–70 rpm	0–25°	Rock/Rotate, Timed Auto-Switch	4.5kg
83041340	Rotator Genie™ (SI-2106AU)	Analog	5–35 rpm	—	—	Rotate	2.25kg
83041381	Roto-Shake Genie™ RT22 (SI-1106AU)	Analog	5–35 rpm	10–70 rpm	10°	Rock, Rotate, Shake with Accessories	4.5kg



MULTI-FUNCTION INCUBATOR GENIE™

Code	Model	Display	Temp. Range	Rotating Speed	Rocking Speed	Stirring Speed	Capacity
83041351	Incubator Genie™ (SI-1406AU)	LED	28°C–75°C	3–35 rpm	6–70 rpm	—	Rotating/Rocking 4.5kg
83041360	Enviro-Genie™ (SI-1206AU)	LCD	4°C–75°C	1–35 rpm	2–70 rpm	1–1000 rpm	Rotating/Rocking 0.5kg. Stirring 2x1L



☎ 02 9575 7512
✉ sales@capellascience.com.au
🌐 capellascience.com.au

**CONTACT CAPELLA SCIENCE ON 02 9575 7512 OR OR
SALES@CAPELLASCIENCE.COM.AU TO BUY NOW**

Freight & GST applies. Valid for purchase orders prior to 30.09.26
Product specifications & availability may be subject to change without prior notice.